

FEATHER® DISPOSABLE SKIN GRAFTING BLADE

Feather disposable skin grafting blade is made of stainless steel. Its sharp edge offers smooth shaving of burns, and prevents snapping.

Skin grafting blade is for use with the following knives:
ROSENBERG, BRAITHWAITE

Dimension (mm) : 157 (L) X 18(W) X 0.229(T)

Material : Stainless Steel

Packaging : Individual packed in aluminum foil and sterilized by Gamma radiation at 25 KG.

Supplied in box of 20 pcs



JAPAN

