



BRANDNEW DIALYSIS MACHINE

BRANDNEW SWS DIALYSIS MACHINE

SWS medical group is focused on extracorporeal blood treatment for 19 years and has become the only company in the extracorporeal blood treatment industry in China with the whole industry chain through technological innovation and continuous improvement.

With association and technology transfer from SWS Medical, we have started assembling SWS Dialysis Machines in India.

Hemodialysis, Off-line Hemofiltration (Off-line HF), Off-line Hemodiafiltration (Off-line HDF), Isolated Ultrafiltration, Sequential Ultrafiltration, Hemoperfusion, etc.

Technical Features:

- 1. Multifunction, Individual Therapy.**
- 2. The Advanced Volume Balance Control System Enables Accurate Control of Ultrafiltration.**
- 3. The Diversified Treatment Modes Can Satisfy Various Clinical Needs.**
- 4. Reliability And Safety Are Ensured By A Dual-Control-Supervising System.**

KEY FEATURES:

- 12.1-inch LCD touch screen
- Kt/v calculation
- Sodium and UF profiles
- HD, IUF, Sequential dialysis, HP
- Venous level detector
- BPM
- Bicarbonate cartridge / Bibag
- Endotoxin Filter (open system)

DIFFERENT THAN OTHERS:

- Battery backup 40 min – 24V Battery
- Data connection software
- Leakage sensors in Hydraulics
- IPX1 standard - Dripping water
- Shall have no harmful effect
- Power Consumption = < 2KVA

- Two years warranty
- Inbuilt Multiple Language Option

COUNTRY OF ORIGIN FOR MAJOR PARTS OF OUR MACHINE

- Valves - made in Germany
- Gear pumps - made in Italy
- Electronic boards chips - made in Japan, assembled in China
- Pressure sensors - made in USA
- Blood pressure monitor - made in USA
- Conductivity sensor - made in Germany
- Heater - made in Italy
- Blood pump - made in Germany
- Pipes - made in Japan
- Silicone tubes (Dow Corning) - made in USA

